



Universal Hip™

Product Mentorship Program

Wet Lab Evaluation Pathway



Our History

BioMedtrix was founded in 1989 with the objective of designing, developing, and manufacturing state-of-the-art veterinary orthopedic implants. Our continuing mission is to provide quality in all aspects of product development, manufacturing, and customer service. Through research sponsorships and collaborations with the world's foremost surgeons, BioMedtrix continues to support the development of new programs to address veterinary needs.

BioMedtrix is represented in over 95% of the US veterinary teaching universities, selected international universities, and specialty clinics throughout the world. We have over 400 clients, some of which have been using BioMedtrix products for over 35 years. This level of support from our customers provides us the opportunity to continue product development.

Research and Development

With over 36 years of development experience in human and veterinary implants, BioMedtrix has developed proprietary processes, and has patents on several designs addressing clinical needs in the canine hip, knee, elbow and trauma. Our research and testing approaches have demonstrated improved product development cycles with designs that promote long term patient mobility.

Quality

BioMedtrix operational policies maintain quality standards comparable to those used by human orthopedic manufacturers. We adhere to all ASTM specifications for implant materials and maintain standards for articulating surfaces of implants which meet or exceed those specified for human implants.

Advisory Program

BioMedtrix maintains relationships with key veterinarians who can provide product development guidance, clinical feedback on the use of these implants, and detailed answers to technical questions. These surgeons are specialists in veterinary orthopedics and cumulatively have over 300 years' experience in veterinary medicine.

Movora

Unleashing the Best Veterinary Outcomes. Together. Through our family of brands, Movora is the premier provider of veterinary orthopedic-related products, education, and innovation worldwide. From infection control through science-driven implant design, to precision instruments and surgical options, we are focused on providing products that help ensure the best veterinary outcomes. BioMedtrix products can be purchased online through movora.com.

Universal Hip

Product Mentorship Program

Wet Lab Evaluation Pathway

Table of Contents

Program Overview	1
BioMedtrix Product Mentorship Program	2
Wet-Lab Evaluation Pathway	4
Process	4-6
Wet Lab Case Evaluation Pathway Completion	6
Fees	7

Program Overview

BioMedtrix Product Mentorship Program

Completion of the BioMedtrix Universal Hip Product Mentorship Program is required before surgeons can purchase implants for clinical use. There are two mentorship programs offered: one that provides approval to purchase implants for standard and large breed dogs, and another that provides approval to purchase implants for small breed dogs and cats.

Each Mentorship Program consists of two parts: a Case Evaluation followed by Clinical Case Mentoring. Within Case Evaluation, surgeons choose one of three pathway options to complete the program. This document covers the Wet Lab Evaluation Pathway, while separate documents are available for the Mentor Evaluation Pathway and the Surgeon Assist Evaluation Pathway.

Surgeons must complete the Universal Hip Standard and Large Breeds Mentorship Program through the Case Evaluation before beginning the Universal Hip Small Breeds Mentorship Program.

There are two specialized mentorship programs offered:

Universal Hip Standard and Large Breeds Mentorship Program	Universal Hip Small Breeds Mentorship Program
<p>Designed for the use of the BioMedtrix Universal Hip System in standard and large breed dogs.</p> <p>Implant Sizes: BFX® Cups 18 - 42 mm CFX® Cups 18 - 31 mm BFX Stems #4 - 13 CFX Stems #4 - 10</p>	<p>Designed for the use of the BioMedtrix Universal Hip system in small breed dogs and cats.</p> <p>Implant Sizes: BFX Cups 10 -18 mm CFX Cups 10 - 18 mm BFX Stems #1.0 - 4 CFX Stems #1 - 4</p>

These programs ensure surgeons receive the necessary training and expertise to perform procedures with confidence and precision.

Both Mentorship Programs follow the same structure and consist of two parts:

Part 1	Case Evaluation	Part 2	Clinical Case Mentoring
---------------	-----------------	---------------	-------------------------

Participation in either Mentorship Program must commence within 12 months of attending the corresponding workshop. Case submissions beyond the 12-month period will require re-attendance of the respective workshop.


Case Evaluation

The Case Evaluation stage of the Mentorship Program offers three potential pathways to practice the procedure, gain confidence, and validate a surgeon's comprehension of the concepts required for successful clinical results. One of the following three pathways may be selected:

Wet Lab Evaluation Pathway	Surgeon Assist Evaluation Pathway	Mentor Evaluation Pathway
<p>Remote Surgeons perform Wet Lab cases in their own clinic.</p> <p>———— or ————</p> <p>In-Person Surgeons perform Wet Lab cases by attending an In-Person Mentorship Opportunity event offered in various locations several times per year.</p> <p>The Wet Lab Evaluation Pathway involves performing Total Hip Replacement surgery on cadavers using the Universal Hip System and submitting pre- and post-operative radiographs, along with a completed Case Evaluation Form for assessment by a BioMedtrix Clinical Advisor.</p> <p><i>Complete pathway details are included within this document.</i></p>	<p>The Surgeon Assist Evaluation Pathway offers the opportunity to continue learning Total Hip Replacement with a BioMedtrix Instructor at the participating surgeon's own veterinary hospital.</p> <p>The instructor will visit the participating surgeon's practice to provide personal coaching and guidance for clinical cases. A maximum of two (2) clinical cases should be performed per day to optimize the learning experience over a negotiated number of days. Note: It may be necessary for the instructor to return to provide guidance for additional cases in order to successfully complete this pathway.</p> <p><i>This pathway is outlined in the Universal Hip Product Mentorship Program: Surgeon Assist Evaluation Pathway document.</i></p>	<p>The Mentor Evaluation Pathway is intended for surgeons who have been assisting an experienced THR surgeon who is proficient in the procedure and recognized by BioMedtrix. Approval of a proposed mentor requires submission of the surgeon's name to biomedtrixmentorship@movora.com. Once approved, the mentor must provide a signed agreement confirming their willingness to support the new surgeon during training. Upon completion of pathway requirements, the mentor will provide a letter acknowledging the surgeon's experience and capability to perform the procedure independently, with consistently successful results.</p> <p><i>This pathway is outlined in the Universal Hip Product Mentorship Program: Mentor Evaluation Pathway document.</i></p>

Program enrollment, pathway selection, and case submission are managed through the Movora Education Platform. Cases are evaluated anonymously by a BioMedtrix Clinical Advisor and must include appropriate pre- and post-operative radiographs, accompanied by a completed Case Evaluation Form. All program details are available on the Movora Education Platform; however, this document provides specific information regarding the Wet Lab Evaluation Pathway.

BioMedtrix Universal Hip™ Education Pathway



Completion of the Educational Pathway is required to purchase implants for clinical cases.

* A limited number of instrument sets are available on a rental basis, subject to availability. Instrument sets may also be purchased.

** Three pathways to completion are available. Delegates select their preference during program enrollment.

*** Following completion of the Case Evaluation step of the Mentorship Program (3a and 6a), participating Surgeons are eligible to purchase a limited number of implants for their first two Clinical Cases (3b and 6b).

➔ After successfully completing the Case Evaluation step of the Mentorship Program (3a), surgeons have the option to directly continue with the BioMedtrix Universal Hip Small Breeds Total Hip Replacement Online Course.

Wet-Lab Evaluation Pathway

The Wet Lab Evaluation Pathway involves performing Total Hip Replacement surgery on canine cadavers using the Universal Hip System and submitting pre- and post-operative radiographs for assessment by a BioMedtrix Clinical Advisor. These practice procedures may be performed at a Movora sponsored “In-Person Mentorship Opportunity,” or at the delegate’s own clinic.

Process

Step 1

Enroll in a Wet Lab Mentorship Program on the Movora Education Platform

https://education.movora.com/visitor_class_catalog

Remote Option	In-Person Option
<p>A</p> <p>Wet Lab Evaluation Pathway: BioMedtrix Universal Hip Mentorship Program – Standard and Large Breeds</p> <p>———— or ————</p> <p>B</p> <p>Wet Lab Evaluation Pathway: BioMedtrix Universal Hip Mentorship Program – Small Breeds</p>	<p>A</p> <p>In-Person Mentorship Opportunity: BioMedtrix Universal Hip Standard and Large Breeds</p> <p>———— or ————</p> <p>B</p> <p>In-Person Mentorship Opportunity: BioMedtrix Universal Hip Small Breeds</p>

Following Registration, you are automatically enrolled in both stages of the Mentorship Program: **Case Evaluation** and **Clinical Case Mentoring**.

Case Evaluation:

Step 2

Acquire a Universal Hip Instrument Set through purchase or the Loaner Program, along with a set of Practice Implants.

- Wet Lab Loaner Instrumentation and Implants may be requested through the Education Platform once a pathway is selected.

- To receive a quote to purchase contact your dedicated Orthopedic Sales Consultant
▶ <https://go.movora.com/Representatives>



Step 3

Complete a minimum number of cases based on the Mentorship Program Selected

Universal Hip Standard and Large Breeds Mentorship Program

Four (4) cadaver cases including three (3) BFX®, and one (1) CFX®, or Hybrid case using the practice implants. If submitting a hybrid case, the stem must be CFX.

Universal Hip Small Breeds Mentorship Program

Two (2) cadaver cases including one (1) BFX, and one (1) CFX, or Hybrid case using the practice implants. If submitting a hybrid case, the stem must be CFX.

Step 4

Submit radiographs and a completed Case Evaluation Form for each case according to the criteria outlined in the THR Case Submission document. Each case will be evaluated and designated “Approved” or “Not Approved” with accompanying explanations.

Note: No more than two (2) cadaver cases may be completed or submitted for evaluation at one time. This limitation ensures the evaluator can thoroughly review and critique each submission while providing recommendations to address any concerns. Identified concerns should then be addressed in submissions of the next two (2) cadaver cases. If more than two (2) cases are submitted simultaneously, only the first two (2) will be evaluated. Participating surgeons must await the completed evaluation of the first two (2) cases before proceeding with additional submissions.

Upon successful completion of the Wet Lab Case Evaluation criteria, the participating surgeon becomes eligible to purchase implants for their initial two (2) clinical cases. As part of the second phase of the Mentorship Program — Clinical Case Mentoring — surgeons are again required to submit radiographs along with a completed Case Evaluation Form.

Clinical Case Mentoring:

Step 5

Proceed with clinical cases and submit radiographs along with a completed Case Evaluation Form for the first two solo clinical cases. Each case will be evaluated and designated “Approved” or “Not Approved” with accompanying explanations.

Wet Lab Case Evaluation Pathway Completion

Successful completion of the Wet Lab Case Training Pathway and the ability to purchase implants for clinical cases occurs when a surgeon achieves each of the following:

Universal Hip Standard and Large Breeds Mentorship Program	Universal Hip Small Breeds Mentorship Program
<ul style="list-style-type: none"> • Receives one (1) Approved CFX or Hybrid wet lab evaluation case. If submitting a hybrid case, the stem must be CFX. • Receives one (1) or more Approved BFX Wet Lab evaluation • Receives an evaluation designation of Approved for the final case submitted for each the CFX/Hybrid and BFX, regardless of prior submitted approved cases, thus demonstrating the ability to correct previously identified issues 	<ul style="list-style-type: none"> • Receives one (1) Approved CFX or Hybrid (CFX Stem/BFX Cup) wet lab evaluation • Receives one (1) Approved BFX Wet Lab evaluation • Receives an evaluation designation of Approved for the final case submitted for each the CFX and BFX, regardless of prior submitted approved cases, thus demonstrating the ability to correct previously identified issues

Surgeons who successfully complete the Mentorship Program—by fulfilling both the Wet Lab Case Evaluation and Clinical Case Mentoring requirements—will be awarded a Universal Hip Product Mentorship Program Certificate of Completion, issued by BioMedtrix. They then have ability to purchase implants for clinical cases.

After completion of the Mentorship Program, Clinical Consultations are available upon request for a fee.

Fees

Wet Lab Evaluation

Remote:

Universal Hip Standard and Large Breeds Mentorship Program – 1500 USD

- Includes the evaluation of 6 cases, which can be a combination of wet lab and clinical cases.
- The Wet Lab Case Evaluation must be completed before beginning Clinical Case Mentoring.
- Any additional cases needed to complete the Mentorship Program will be invoiced at 250 USD per case.

Equipment Fees:

- Universal Hip Instrument Loaner **1200 USD**
- Universal Hip Practice Implants **500 USD**
- Applicable shipping fees will be invoiced

Universal Hip Small Breeds Mentorship Program – 1000 USD

- The Wet Lab Case Evaluation must be completed before beginning Clinical Case Mentoring.
- Any additional cases needed to complete the Mentorship Program will be invoiced at 250 USD per case.

Equipment Fees:

- Universal Hip Small Breeds Instrument Loaner **1200 USD**
- Universal Hip Small Breeds Practice Implants **500 USD**
- Applicable shipping fees will be invoiced

Loaner implants and instruments are to be returned following the guidelines outlined in the **Loaner Return** document.

In Person:

Advertised price includes all supplies, case evaluation, and clinical case mentoring.

See Movora Education Platform for current rates.

For more information regarding the Product Mentorship Program, please email biomedtrixmentorship@movora.com.

Outside the United States your representative can inform you of program fees.



Movora
Local contact: www.movora.com/contact-us

Subject to technical changes.
©Movora 2025. All rights reserved.
Doc. Nr.: MMKT-ED-0066-A